
COVID-19 Antigen Detection Kit

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
Dr. M.J. van de Velde

CE Certification – Declaration of Conformity

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DECLARATION OF CONFORMITY

Regarding In Vitro Diagnostic Directive (98/79/EC)

Manufacturer: New Gene (Hangzhou) Bioengineering Co., Ltd.

Address: Room 1606, Floor 16, Building 5, 688 Bin'an Road, Changhe Street, Binjiang District, Hangzhou City, Zhejiang Province, P. R. China

EC Representative: SUNGO Europe B.V.

Address: Olympisch Stadion 24, 1076DE Amsterdam, Netherlands

Product Name: COVID-19 Antigen Detection Kit

Specification: 25Tests/Box 1Test/Box

Classification: Others (IVDD)

Conformity Assessment

Procedure: Annex III of In Vitro Diagnostic Directive (98/79/EC)

We herewith declare that the above-mentioned products meet the requirements of In Vitro Diagnostic Directive (98/79/EC) and the following harmonized standards.

EN 23640:2015

EN 13640:2002

EN 980:2016

EN 13641:2002

EN ISO 14971:2019

EN ISO 18113-1 2011

EN 13612:2002

EN ISO 18113-4 2011

Signature:

Name/ Position: Mingfu Li / General Manager

On behalf of SUNGO Europe office, I confirmed we are EU REP of the company who issue this document.

Date: 29/09/2020

Authorized Signature (S)

Place: Hangzhou, Zhejiang, China

By Royal Charter

Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that: New Gene (Hangzhou)
Bioengineering Co., Ltd.
Room 1606, 16th Floor, No.5 Building
688 Bin'an Road
Binjiang District
Hangzhou
Zhejiang
310052
China

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Holds Certificate No: **MD 729179**

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

Design and Development, Manufacture and Distribution of In-vitro Diagnostic Rapid Test Kit of Drug Abuse, Manufacture and Distribution of In-vitro Diagnostic Rapid Test Kit of Infectious Diseases.

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For and on behalf of BSI:

Gary E Slack, Senior Vice President - Medical Devices

Original Registration Date: 2020-07-27

Latest Revision Date: 2020-07-27

Effective Date: 2020-07-27

Expiry Date: 2023-07-26

Page: 1 of 1

...making excellence a habit.™

Included in the “Export Allow List” 出口白名单企业

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91330108MA2H3RX57K

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New Gene (Hangzhou) Bioengineering Co., Ltd.

Certification for Safe Transport 货物运输鉴定

NO.2020107244

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Certification
for Safe Transport of Chemical Goods

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Sample Name: COVID-19 Antigen Detection Kit

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Shanghai Research Institute of Chemical Industry Testing Co., Ltd

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Ministry of Health, Welfare and Sport
CIBG
P.O. Box 16114
2500 BC The Hague
THE NETHERLANDS

STATEMENT

The undersigned herewith declares that according to the Decree on In-Vitro Diagnostics, which is based on the European Directive 98/79/EC concerning in-vitro diagnostic medical devices,

SUNGO Europe B.V.
Olympisch Stadion 24
1076 DE Amsterdam
THE NETHERLANDS

acts as authorised representative of the manufacturer.

The manufacturer:

New Gene (Hangzhou) Bioengineering Co., Ltd.
Room 1606, Floor 16, Building 5, 688 Bin'an Road, Changhe Street, Binjiang District,
Hangzhou City, Zhejiang Province
CHINA

is authorised to manufacture and/or supply the medical device/devices mentioned below:

COVID-19 / Influenza A / Influenza B Detection Kit
COVID-19 Antibody / Antigen Detection Kit
COVID-19 Antigen Detection Kit
COVID-19 Neutralizing Antibody Detection Kit
Novel Coronavirus Ribonucleic Acid Detection Kit

This device/these devices may be placed on the Dutch market and on the markets of the other Member States of the European Union, and be exported to non-EU Member States. This free sale certificate may only be used for export outside the European Union.

The present statement is drawn up at the request of the interested party in order to be submitted to the Health Authorities of **INDIA**.

This statement is valid until May 26, 2022.

The Hague, October 20, 2020

On behalf of the Minister for Medical Care and Sport
Farmatec | CIBG

Dr. M.J. van de Velde

Mr. M.J. van de Velde
Head of Department

Our reference: 20204982
Certificate number: 29432

COVID-19 Antigen Detection Kit

Package Insert

Cat: COVID-19-NG08 **Specimens:** Sputum
Version: 03 **Effective Date:** 2020-10

For professional and in vitro diagnostic use only.

PRODUCT NAME

COVID-19 Antigen Detection Kit

PACKING

1 piece/bag, 25 pieces/box.

INTENDED USE

This product is suitable for the qualitative detection of novel coronavirus in sputum sample. It provides an aid in the diagnosis of infection with novel coronavirus.

SUMMARY

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

PRINCIPLE

The COVID-19 Antigen Detection Kit is an immunochromatographic membrane assay that uses highly sensitive monoclonal antibodies to detect nucleocapsid protein from SARS-CoV-2 in throat swab. The test strip is composed of the following parts: namely sample pad, reagent pad, reaction membrane, and absorbing pad. The reagent pad contains the colloidal-gold conjugated with the monoclonal antibodies against the nucleocapsid protein of SARS-CoV-2; the reaction membrane contains the secondary antibodies for nucleocapsid protein of SARS-CoV-2. The whole strip is fixed inside a plastic device. When the sample is added into the sample well, conjugates dried in the reagent pad are dissolved and migrate along with the sample. If SARS-CoV-2 antigen presents in the sample, a complex formed between the anti-SARS-2 conjugate and the virus will be captured by the specific anti-SARS-2 monoclonal antibodies coated on the test line region (T). Absence of the T line suggests a negative result. To serve as a procedural control a red line will always appear in the control line region (C).

indicating that proper volume of sample has been added and membrane wicking has occurred.

COMPOSITION


1. Test Card
2. Sample Extraction Tube
3. Cotton Swab
4. Paper Cup

STORAGE AND STABILITY


1. Store as packaged in the hermetic bag at the temperature (2-30°C or 38-86°F) and avoid direct sunshine. The kit is stable within the expiration date printed on the labeling.
2. Once open the hermetic bag, the test should be used within one hour. Prolonged exposure to hot and humid environment will cause product deterioration.
3. The lot number and the expiration date are printed on the labeling.

TEST PROCEDURE

Allow the test device and specimens to restore to room temperature (15-30°C or 59-86°F) prior to testing.

Sampling Steps

Detection Steps

1. Use the cotton swab to pick up 10-50 mg sputum samples (equivalent to the size of a match head). Open the cap of sample extraction tube, break the swab tip into the tube. Close the sample extraction tube and shake to mix the sample completely. Leave the swab in the extraction tube for one minute.
2. Take the test card from the packaging bag, place it on a table, cut off the

protrusion of the collection tube, and add 3 drops of the sample solution into the sample loading hole vertically.


3. Read the result in 15 minutes. The result is considered inaccurate and invalid after 30 minutes.

INTERPRETATION OF RESULTS

Positive(+): Both of T and C lines are appeared in 15minutes.

Negative(-): C line is appeared while no T line appeared in 15 minutes after the sample is loaded.

Invalid: As long as the C line does not appear, it indicates that the test result is invalid, and should retest the sample with another test card.

NOTES

1. The COVID-19 Antigen Detection Kit is applicable to sputum samples. Blood, serum, plasma, urine and other samples may cause abnormal results.
2. Please make sure that a proper amount of sample is added for testing. Too much or too little sample amount may cause deviations in results.
3. For positive judgement, it can be confirmed as soon as both T and C line appeared. That may take 3-15 minutes after the sample is loading. For negative judgement, please wait for 15 minutes after sample loading. The result is invalid after 30 minutes after sample loading.
4. This product is disposable. DO NOT recycle.
5. Dispose of used products, samples, and other consumables as medical wastes under relevant regulations.
6. If the test line or control line is out of the test window, do not use the test card. The test result is invalid and retest the sample with another one.

INDEX OF SYMBOLS

	Consult instructions for use		Tests per kit		Authorized Representative
	For in vitro diagnostic use only		Use by		Do not reuse
	Store between 2-30°C		Lot Number		Catalogue number

New Gene (Hangzhou) Bioengineering Co., Ltd.
Room 1606, Floor 16, Building 5, 688 Bin'an Road,
Changhe Street, Binjiang District, Hangzhou City,
Zhejiang Province, P. R. China

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COVID-19 Antigen Detection Kit

Package Insert

Cat: COVID-19-NG08 **Specimens:** Throat Swab
Version: 02 **Effective Date:** 2020-10

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PRODUCT NAME
COVID-19 Antigen Detection Kit

PACKING
1 piece/bag, 25 pieces/box.

INTENDED USE
This product is suitable for the qualitative detection of novel coronavirus in throat swab sample. It provides an aid in the diagnosis of infection with novel coronavirus.

SUMMARY
The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

PRINCIPLE
The COVID-19 Antigen Detection Kit is an immunochromatographic membrane assay that uses highly sensitive monoclonal antibodies to detect nucleocapsid protein from SARS-CoV-2 in throat swab. The test strip is composed of the following parts: namely sample pad, reagent pad, reaction membrane, and absorbing pad. The reagent pad contains the colloidal-gold conjugated with the monoclonal antibodies against the nucleocapsid protein of SARS-CoV-2; the reaction membrane contains the secondary antibodies for nucleocapsid protein of SARS-CoV-2. The whole strip is fixed inside a plastic device. When the sample is added into the sample well, conjugates dried in the reagent pad are dissolved and migrate along with the sample. If SARS-CoV-2 antigen presents in the sample, a complex formed between the anti-SARS-2 conjugate and the virus will be captured by the specific anti-SARS-2 monoclonal antibodies coated on the test line region (T). Absence of the T line suggests a negative result. To serve as a procedural control a red line will always appear in the control line region (C) indicating that proper volume of sample has been added and membrane wicking has occurred.

COMPOSITION


1. Test Card
2. Sample Extraction Tube
3. Throat Swab

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SAMPLE COLLECTION


1. Use the throat swab provided in the kit to swab over the lateral and posterior walls of pharynx, as well as the intratonsillar cleft.
2. Open the cap of sample extraction tube, break the swab tip into the tube. Close the sample extraction tube and shake to mix the sample completely. Leave the swab in the extraction tube for one minute.

TEST PROCEDURE


Allow the test device and specimens to equilibrate to temperature (15-30°C or 59-86°F) prior to testing.

1. Take the test card from the packaging bag, place it on a table, cut off the protrusion of the collection tube, and add 3 drops of the sample solution into the sample loading hole vertically.
2. Read the result in 15 minutes. The result is considered inaccurate and invalid after 30 minutes.

INTERPRETATION OF RESULTS

Positive(+): Both of T and C lines are appeared in 15 minutes.
Negative(-): C line is appeared while no T line appeared in 15 minutes after the sample is loaded.
Invalid: As long as the C line does not appear, it indicates that the test result is invalid, and should retest the sample with another test card.

NOTES

1. The COVID-19 Antigen Detection Kit is applicable to throat swab sample. Blood, serum, plasma, urine and other samples may cause abnormal results.
2. Please make sure that a proper amount of sample is added for testing. Too much or too little sample amount may cause deviations in results.
3. For positive judgement, it can be confirmed as soon as both T and C line appeared. That may take 3-15 minutes after the sample is loading. For negative judgement, please wait for 15 minutes after sample loading. The result is invalid after 30 minutes after sample loading.
4. This product is disposable. DO NOT recycle.
5. Dispose of used products, samples, and other consumables as medical wastes under relevant regulations.
6. If the test line or control line is out of the test window, do not use the test card. The test result is invalid and retest the sample with another one.

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FAQ NEWGENE Coronavirus Antigen Test

1 What is Coronavirus Antigen test?

Coronavirus Antigen test detects the virus itself by the specific proteins (i.e. Antigens) on its surface, so as to identify the existence of the virus.

NEWGENE technician is preparing the solution

2 What is the difference between Coronavirus Antigen test and Antibody test?

Antibody test is to detect the molecules (i.e. antibodies) that people produce after being infected with the virus. Antibody may take several days to produce after infection, and usually remain in the blood for several weeks after the patient recovers. That is, the cured patient can also be tested positive for a longer period of time. Therefore, the antibody test is useful in the diagnosis of COVID-19, but with limitations.

The Antigen test is to detect the virus, which directly relates to the infection itself.

3 What is the difference between Antigen test and nucleic acid test?

Antigen test is based on immune reaction, while nucleic acid test is based on genetic information. Therefore, Antigen test can be much faster than nucleic acid test (PCR) (15 minutes vs 4 hours). Certainly, the cost of speed is the sensitivity.

That is, Antigen test is not as sensitive as nucleic acid test. PCR can detect micro amount of SARS-CoV-2 virus due to the amplification process of virus genetic information. This is also why PCR test often takes several hours and may cause aerosol contamination.

4 Advantages of Antigen test and precautions in application

Except for advantages of rapid test, it does not require laboratory processing and medical professionals to operate. Therefore, it is very suitable for large-scale general screening when lacking laboratory conditions and medical professionals. However, if human body contains a very low amount of virus, and with no obvious clinical symptoms, they are likely to be tested negative. The ability of such infected people to spread the virus is also much lower than that of ordinary patients. Therefore, for asymptomatic infected people, it is recommended to adopt nucleic acid for detection.

Antigen test helps to quickly identify people with high levels of infection (those who are most likely to infect others). Therefore, Antigen test can help prevent the spread of the pandemic, by identifying those who are most likely to spread the disease and isolating them from others. Thus, shift the focus to identify the most infectious people.

NEWGENE technicians are making samples


5 Precautions for Antigen test operations

Antigen test is to detect the virus itself, so sample selection is very important. We recommend sputum sample, since sputum is the secretion from respiratory tract where viruses cluster the most, so sputum is relatively easy for detection.

Because of the uneven distribution of the virus on the nasopharyngeal mucosa appearance, it is very likely that no virus or very small amount of virus is taken while sampling and lead to false negative result, when sampling via throat swab. Therefore, the accuracy will be lower than that of sputum.

In addition, studies have shown that after people are infected, the viruses can be detected in the digestive tract, especially in the intestinal stool. Therefore, if convenient (e.g. in hospital or at home), stool sample is another option for testing; Saliva samples mixed with sputum also have a certain probability of detection, although the amount of virus contained is lower than pure sputum.

For false negative testing result (PCR positive, Antigen test negative), if the patient has asymptomatic infection, it is recommended to take nucleic acid as standard, and if patients cough or have other symptoms, it's suggested to re-sampling and test again. Besides, S protein is selected as the testing object in the NEWGENE's Antigen test. S protein shows good specificity and usually does not cause false positive results.

NEWGENE technicians are testing different concentrations of standards

6 Others

The WHO has always been a supporter of Antigen testing. Experts pointed out that compared with conducting an accurate test every two weeks, conducting a relatively insensitive quick test twice a week can more effectively contain the spread of SARS-CoV-2. That is, the focus of the virus detection should be on identifying those who are likely to spread SARS-CoV-2, rather than locating anyone who is infected with SARS-CoV-2.

Delhi is the first place in India to start using rapid Antigen testing in June. By mid-July, the number of cases there had been reduced, and the daily death toll had also stabilized, indicating that Antigen testing played a certain role in controlling the spread of the virus. Although after cases declined, India began to lift the blockade restrictions and cause the infection to rise again.

The Philippine Association of Microbiology and Infectious Diseases has issued interim guidelines for clinicians and medical staff, stating that Antigen testing can be used as a substitute for PCR to diagnose coronavirus infection in people with symptoms in the early stage, especially during the acute infection period.

Delhi is the first place in India to start using rapid Antigen testing in June. By mid-July, the number of cases there had been reduced, and the daily death toll had also stabilized, indicating that Antigen testing played a certain role in controlling the spread of the virus. Although after cases declined, India began to lift the blockade restrictions and cause the infection to rise again.

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