COVID-19 Antigen Rapid Test Clinical Sensitivity and Specificity Study Report

1. Objective

The CLUNGENE® COVID-19 Antigen Rapid Test (hereinafter referred to as the CLUNGENE Device) manufactured by Hangzhou Clongene Biotech Co., Ltd. is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasopharyngeal swab and oropharyngeal swab from individuals who are suspected of COVID-19 by their healthcare provider.

This study is intended to evaluate the clinical performance, between the CLUNGENE Device and the comparator RT-PCR assay.

2. Method

A study of 285 direct nasopharyngeal swabs was performed. The specimens were collected from symptomatic patients suspected of COVID-19 at 3 locations and tested at a single central laboratory.

Two nasopharyngeal swabs were collected from individual symptomatic patients (within 7 days of onset) who were suspected of COVID-19. At all locations, one nasopharyngeal swab was immediately frozen at -70°C for later testing, and the other swab was eluted in 3 mL viral transport media and tested with RT-PCR assay for detection of SARS-CoV-2. A total of 285 retrospective nasopharyngeal swab specimens within a pre-specified date range were selected and then tested by the CLUNGENE Device in a blinded fashion. The operators were blinded to the RT-PCR test results.

The positive percent agreement (PPA) was calculated as $100\% \times (\text{True Positive} / [\text{True Positive} + \text{False Negative}])$. The negative percent agreement (NPA) was calculated as $100\% \times (\text{True Negative} / [\text{True Negative} + \text{False Positive}])$. The 95% (two-sided) confidence interval (CI) was calculated using the Wilson Score Method.

3. Comparator method

Real-Time Fluorescent RT-PCR Kit for Detecting SARS-CoV-2, manufactured by BGI Genomics Co. Ltd. This product has got CE, NMPA certification and FDA Emergency Use Authorized. A specimen is positive for SARS-CoV-2 if the Ct value of ORF1ab gene is not higher than 37.

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4. Enrollment criteria (inclusion/exclusion criteria)

4.1 Inclusion criteria

- Individual symptomatic patients (within 7 days of onset) who were suspected of COVID-19.

4.2 Exclusion criteria

- Unable to obtain samples of information needed for the experiment
- Samples that have been contaminated or contaminated during sample storage
- Samples with inappropriate storage conditions

5. Result

The results are summarized in the following table.

COVID-19 Antigen		RT-PCR		Total
		Positive	Negative	
CLUNGENE®	Positive	64	0	64
	Negative	6*	215	221
Total		70	215	285

Positive Percent Agreement (PPA)= 91.4% (64/70), (95%CI: 82.5%~96.0%)

Negative Percent Agreement (NPA) = 100% (215/215), (95%CI: 98.2%~100%)

* The 6 discordant specimens (CLUNGENE Device Negative/ Comparator RT-PCR assay Positive) had Ct Values of 34, 36, 35.5, 34, 35 and 33.

The PPA is 98.5% (64/65) (95%CI: 91.8%~99.7%) with specimens of a Ct count ≤ 33 .

6. Conclusion

Taken together, the CLUNGENE Antigen Rapid Test had a positive percent agreement (sensitivity) of 91.4% (95% CI: 82.5%~96.0%) and negative percent agreement (specificity) of 100% (95% CI: 98.2%~100%). The PPA is 98.5% (64/65) (95%CI: 91.8%~99.7%) with specimens of a Ct count ≤ 33 .

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