

Information on Virus Mutations

Recently, a new variant of the SARS-CoV-2 virus, the causative agent of COVID-19, has become known. This variant, named Omicron (B.1.1.529), was first found in South Africa and shows a rapid spread and high number of mutations [1]. This raises the question of what effect the mutations have on detection by the Testsealabs Antigen Rapid Test.

As things stand, detection of COVID-19 disease by the Testsealabs Antigen Rapid Test is reliable even with the Omicron variant. Even though the variant shows a large number of mutations, these mutations are located in the part called spike protein on the surface of the virus [2]. This part helps the virus to enter cells, but at the same time makes it recognizable to the immune system.

Detection by the Testsealabs Antigen Rapid Test is based on another protein called the nucleocapsid protein. This shows only a few mutations in the Omicron variant, some of which also occur in other variants such as Alpha, Gamma or Lambda. The mutations are also located outside the crucial components of the protein, which are a prerequisite for recognition. The antibodies used in the Testsealabs Antigen Rapid Test recognize regions (epitopes) in the RNA-binding domain in the range of amino acids 47-173. This part of the protein is not altered by mutations, so that, according to the current status, reliable detection should be possible. This is also supported by preliminary results with antigen rapid tests from other manufacturers [3].


Harry He, CEO
Hangzhou Testsea Biotechnology Co.,Ltd


[1] Classification of Omicron (B.1.1.529): SARS-CoV-2 Variant of Concern, World Health Organization, 2021
 [2] European Centre for Disease Prevention and Control. Implications of the emergence and spread of the SARS-CoV-2 B.1.1.529 variant of concern (Omicron), for the EU/EEA. 26 November 2021. ECDC: Stockholm; 2021.
 [3] Prof. Dr. Sandra Ciesek, Johann Wolfgang Goethe-Universität Frankfurt, <https://twitter.com/CiesekSandra/status/1464925734741368837?s=20>

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[1] Classification of Omicron (B.1.1.529): SARS-CoV-2 Variant of Concern, World Health Organization, 2021

[2] European Centre for Disease Prevention and Control. Implications of the emergence and spread of the SARS-CoV-2 B.1.1.529 variant of concern (Omicron), for the EU/EEA. 26 November 2021. ECDC: Stockholm; 2021.

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